

# BlueROCKER® Short

BlueROCKER® SHORT is specifically developed to accommodate individuals with a shorter stature, slightly broader shaped tibia, and a wider foot at the metatarsal area. Tibia padding should always be used for patient comfort. This model has a heel height of 7 mm in all sizes.

## Recommended Range Of Application

Foot drop in combination with weak plantar flexors, and in some cases limited knee control. Intended for moderate to severe spasticity and proximal deficits, impairment in multiple leg muscle groups and partial foot amputation. Also for bilateral wearers and balance issues.

## Contraindications

Should not be used when patients present with foot and/or leg ulcers, moderate to severe edema, moderate to severe foot deformities, severe proximal deficits or severe spasticity.

BLUEROCKER® SHORT



Item No.	Product Description	Size	L/R	Color
289951011	BlueROCKER® Short	Small	Left	-
289951012	BlueROCKER® Short	Medium	Left	-
289951013	BlueROCKER® Short	Large	Left	-
289952011	BlueROCKER® Short	Small	Right	-
289952012	BlueROCKER® Short	Medium	Right	-
289952013	BlueROCKER® Short	Large	Right	-
287680011	BlueROCKER® Short, Wrap Around Straps	Small	-	-
287680012	BlueROCKER® Short, Wrap Around Straps	Medium	-	-
287680013	BlueROCKER® Short, Wrap Around Straps	Large	-	-
283910011	SoftKIT™ for ToeOFF® & BlueROCKER® Short	Small	-	Black
283910012	SoftKIT™ for ToeOFF® & BlueROCKER® Short	Medium	-	Black
283910013	SoftKIT™ for ToeOFF® & BlueROCKER® Short	Large	-	Black