

SELECTION® Thumb Orthosis Rigid

Soft and elastic thumb orthosis with extra stability around the thumb structures. The design supports the muscles, keeping the thumb in neutral position and unloads CMC and MPC, but allows an active "pinch" grip. At rest the thumb is kept in neutral position without straining the muscles and tendons. Can be trimmed to the desired shape.

Recommended Range Of Application

Mild arthritis in the first CMC-joint. Protects the CMC-join and the first MPC-joint.

After trauma, in case of hyper mobile thumb and wrist, after arthroplasty surgery.

SELECTION®



Item No.	Product Description	Size	L/R	Color	Circumference
352021010	Selection Thumb Orthosis Rigid	X-Small	Left	Black	< 5in
352021011	Selection Thumb Orthosis Rigid	Small	Left	Black	5in-6in
352021012	Selection Thumb Orthosis Rigid	Medium	Left	Black	6in-6.75in
352021013	Selection Thumb Orthosis Rigid	Large	Left	Black	6.75in-8in
352021014	Selection Thumb Orthosis Rigid	X-Large	Left	Black	> 8in
352022010	Selection Thumb Orthosis Rigid	X-Small	Right	Black	< 5in
352022011	Selection Thumb Orthosis Rigid	Small	Right	Black	5in-6in
352022012	Selection Thumb Orthosis Rigid	Medium	Right	Black	6in-6.75in
352022013	Selection Thumb Orthosis Rigid	Large	Right	Black	6.75in-8in
352022014	Selection Thumb Orthosis Rigid	X-Large	Right	Black	> 8in