

Vission Shoulder Immobilizer

Vission shoulder immobilizer is intended to immobilize and contain the upper limb, elbow, forearm flexed at 90° or 130° as required. Especially suited for shoulder and pectoral girdle lesions. The product is reversible and can be used either on the left or on the right shoulder as required.

Recommended Range Of Application

Shoulder, arm, elbow, and pectoral girdle trauma (fractures, dislocations, and sprains). Shoulder and pectoral girdle surgery. Paralysis of brachial-plexus.



Item No.	Product Description	Size	Color	Length	Chest Circumference
460220010	Vission Shoulder Immobilizer	X-Small	Navy Blue	< 11.75in	26.75"-30.75"
460220011	Vission Shoulder Immobilizer	Small	Navy Blue	< 15in	31.5"-35.5"
460220012	Vission Shoulder Immobilizer	Medium	Navy Blue	< 15in	35.5"-39.25"
460220013	Vission Shoulder Immobilizer	Large	Navy Blue	< 15in	39.25"-44.5"
460220014	Vission Shoulder Immobilizer	X-Large	Navy Blue	< 15in	44.5"-47.25"
460220015	Vission Shoulder Immobilizer	XX-Large	Navy Blue	< 15in	47.25"-55"