

Binky Universal Free Motion Joint

These universal joints, without flexion and/or extension stops, are paired in the production of orthoses from thermoplastics or separately used as a companion free motion joint.

Recommended Range Of Application

Suitable for use by both children and adults in the wrist, elbow, and ankle joint orthoses.

Contraindications

Paralytic syndromes.



Item No.	Length	Width	Thickness	User Weight
696880011	62 mm	12 mm	2, 5 mm	< 35 kg
696880013	76 mm	18 mm	2, 5 mm	< 100 kg