

## SWASH® GO

(Replaces previous SWASH Low Profile. A “Conversion Chart” can be found below under the DOCUMENTS tab.)

Swash is a unique dynamic hip brace that allows controlled variable abduction and dynamic movement in all three planes of motion. The goal of SWASH® is to prevent excessive hip adduction when sitting, standing, and walking as usual in cerebral palsy. However, perhaps the most important benefit the SWASH® offers is guidance for hip alignment.

SWASH GO is the choice for the standing and ambulatory patient, when primary goal is to control scissoring gait and when patient has limited space between the iliac crest and the rib cage. It also encourages greater trunk extension for improved postural positioning.

A complete SWASH orthosis requires 3 parts: a Pelvic Section, Thigh Cuffs, & Uprights. Please order each item separately.

SWASH GO Pelvic section has an aluminum posterior frame that extends from approximately L2/L3 to S1/S2, with an iliac extension.

Uprights are available in two versions: 115° = narrow sitting base / 123° = wide sitting base

### **Recommended Range Of Application**

When primary goal is to control scissoring gait and when patient has limited space between the iliac crest and the rib cage.

### **Measurements**

At the level of the ASIS. Ideal is to fit the orthosis at mid-pelvic girdle to obtain good pelvic stability. See "Product Options" below for Sizing Guide.



Item No.	Product Description	Size	Circumference	Length	Diameter
288530001	GO Pelvic Section	1	40-46,5 cm	-	-
288530002	GO Pelvic Section	2	46-55 cm	-	-
288530003	GO Pelvic Section	3	54-63 cm	-	-
288530004	GO Pelvic Section	4	62-72 cm	-	-
288530005	GO Pelvic Section	5	71-81 cm	-	-
288530006	GO Pelvic Section	6	80-90 cm	-	-
288540801	GO Thigh Cuff	1	21-25 cm	-	-
288540802	GO Thigh Cuff	2	25-29 cm	-	-
288540803	GO Thigh Cuff	3	29-33 cm	-	-
288540804	GO Thigh Cuff	4	33-38 cm	-	-
288540805	GO Thigh Cuff	5	38-44 cm	-	-
288540806	GO Thigh Cuff	6	44-50 cm	-	-
288550801	GO Uprights 115°	1	-	31 cm	8 mm
288550802	GO Uprights 115°	2	-	37.5 cm	8 mm
288550803	GO Uprights 115°	3	-	44 cm	8 mm
288550804	GO Uprights 115°	4	-	48 cm	8 mm
288550805	GO Uprights 115°	5	-	54 cm	8 mm
288550806	GO Uprights 115°	6	-	60 cm	8 mm
288560801	GO Uprights 123°	1	-	31 cm	8 mm
288560802	GO Uprights 123°	2	-	37.5 cm	8 mm
288560803	GO Uprights 123°	3	-	44 cm	8 mm
288560805	GO Uprights 123°	5	-	54 cm	8 mm
288560806	GO Uprights 123°	6	-	60 cm	8 mm