

# SELECTION® Finger Loop

An orthosis that efficiently holds two fingers or more together. A stronger finger assists a weaker one for support, protection and/or increased range of motion. It can also be used as a static progressive splint for an increased range of motion in the DIP joint.

The orthosis is made of 100% durable and quick-drying nylon. Welded with no seams and therefore easy to cut for a personal fit. One or more pieces can be used at the same time.

The semi-prepared Finger Loop sheet, which is 10 cm wide, can easily be trimmed to individual straps in the desired width using scissors.

### Recommended Range Of Application

Support for the fingers after e.g., injury, surgery, stable fractures, metacarpal fractures, chronic pain conditions, sprains, joint wear, and contractures.



Item No.	Product Description	Size	Length	Width	Unit/Pack
352190000	SELECTION® Finger Loop	One Size	15.5 cm	2 cm	5
352200000	SELECTION® Finger Loop Sheet	One Size	15.5 cm	10 cm	1