

SELECTION Radialis Thumb Part

Supporting textile thumb section to be combined with SELECTION® Wrist orthoses 35204, 35205, 35209, 35222, 35223, or 35305. Read the instructions for the main product before use. In case of radial nerve palsy (wrist drop caused by damage to the radial nerve), the thumb section can be combined with the 35217 SELECTION® Radial Finger Part to enable a functional grip while the nerve recovers.



Recommended Range Of Application
For relief and guidance of the thumb, Radial Nerve Palsy.

Measurements
Length of joint between center of MCP1 to joint center of IP1.

Item No.	Product Description	Size	Measurement	L/R	Color
352161011	Radialis Thumb Part	Small	3 - 4 cm	Left	Black
352161013	Radialis Thumb Part	Large	5 - 6 cm	Left	Black
352162011	Radialis Thumb Part	Small	3 - 4 cm	Right	Black
352162013	Radialis Thumb Part	Large	5 - 6 cm	Right	Black