

# SELECTION Radialis Finger Part

Supporting textile finger section to be combined with SELECTION® Wrist orthoses 35204, 35205 or 35209. Read the instructions for the main product before use. In case of radial nerve palsy (wrist drop caused by damage to the radial nerve), the finger section can be combined with the 35216 SELECTION® Radial Thumb Part to enable a functional grip while the nerve recovers.



## Recommended Range Of Application

Radial Nerve Palsy

## Measurements

Circumference MCP 2-5.

Item No.	Product Description	Size	Measurement	L/R	Color
352171011	Radialis Finger Part	Small	17 - 20 cm	Left	Black
352171012	Radialis Finger Part	Medium	20 - 23 cm	Left	Black
352171013	Radialis Finger Part	Large	23 - 27 cm	Left	Black
352172011	Radialis Finger Part	Small	17 - 20 cm	Right	Black
352172012	Radialis Finger Part	Medium	20 - 23 cm	Right	Black
352172013	Radialis Finger Part	Large	23 - 27 cm	Right	Black