

SELECTION Wrist Orthosis Rigid

Selection Wrist Orthosis Rigid has the same features as the soft version but includes a thumb section for increased thumb support and stability. The thumb strap can be individually adjusted. It comes with a pre-shaped plastic spoon on the palmar side and a soft narrow plastic spoon on the dorsal side. This product is intended to be used when the first CMC-joint and wrist needs support and off-loading in cases such as rheumatic disorders, soft tissue injuries or attritional wear.

Made with the airy and comfortable SELECTION material with all the functionality this incorporates.

Recommended Range Of Application

Inflammation, work related strain or when rheumatic and/or degenerative changes are present. Also, for soft tissue injuries and attritional wear.

Contraindications

Non known.

Measurements

Wrist circumference.



Item No.	Product Description	Size	Measurement	L/R
352061010	Rigid	X-Small	< 11 cm	Left
352061011	Rigid	Small	11 - 14 cm	Left
352061012	Rigid	Medium	14 - 17 cm	Left
352061013	Rigid	Large	17 - 20 cm	Left
352061014	Rigid	X-Large	20 - 23 cm	Left
352061015	Rigid	XX-Large	23 - 26 cm	Left
352062010	Rigid	X-Small	< 11 cm	Right
352062011	Rigid	Small	11 - 14 cm	Right
352062012	Rigid	Medium	14 - 17 cm	Right
352062013	Rigid	Large	17 - 20 cm	Right
352062014	Rigid	X-Large	20 - 23 cm	Right
352062015	Rigid	XX-Large	23 - 26 cm	Right