

SELECTION Multiopen

SELECTION MultiOpen is a medium-stability wrist orthosis with a dual locking function. It can be fully opened and then loosely closed with Velcro® and tightened with lacing. The user can choose to open it completely each time the orthosis is donned or doffed, or when only the lacing function is used. It is low cut to follow the hand's natural bend and the area around the thumb allows the thenar muscles to work freely. The two-piece drawstring in combination with the Velcro closure allows for the adjustment of the width for a wider forearm. As the lacing is on the ulnar side, the risk of misalignment of the thumb is reduced. There is also an adjustable thumb band on the volar side. Due to the construction of the material, it is possible to slightly adjust the shape by trimming with scissors for optimal fit.

Made with the airy and comfortable SELECTION material with all the functionality this incorporates. The orthosis is 22 cm long, 4 cm longer than a SELECTION Wrist Soft 35204. It has a removable moldable aluminum spoon on the volar side.

Recommended Range Of Application

Aids in recovery after radius fracture, soft tissue damage, inflammation, wear and tear, work-related overload, and in rheumatic or degenerative changes.

Contraindications

Not to be used for non-healed fractures.

Measurements

Wrist circumference.



Item No.	Product Description	Size	Measurement	L/R
352091010	SELECTION Multiopen	X-Small	<11 cm	Left
352091011	SELECTION Multiopen	Small	11 - 14 cm	Left
352091012	SELECTION Multiopen	Medium	14 - 17 cm	Left
352091013	SELECTION Multiopen	Large	17 - 20 cm	Left
352091014	SELECTION Multiopen	X-Large	20 - 23 cm	Left
352091015	SELECTION Multiopen	XX-Large	23 - 26 cm	Left
352092010	SELECTION Multiopen	X-Small	< 11 cm	Right
352092011	SELECTION Multiopen	Small	11 - 14 cm	Right
352092012	SELECTION Multiopen	Medium	14 - 17 cm	Right
352092013	SELECTION Multiopen	Large	17 - 20 cm	Right
352092014	SELECTION Multiopen	X-Large	20 - 23 cm	Right
352092015	SELECTION Multiopen	XX-Large	23 - 26 cm	Right