

SELECTION Bromsa Wrist Orthosis

Selection Bromsa Wrist Orthosis with Dorsal Stay is intended to restrict movement in the wrist to off-load the hand extensors and support the wrist to reduce pain in conditions such as tennis elbow, hypermobility, rheumatic problems, and unspecific wrist pain. When combined with our Vission EPI Tennis Elbow orthosis, it may help improve grip strength and reduce pain in some patients with tennis elbow. The Selection Bromsa Wrist orthosis has two plastic stays on the dorsal side and an adjustable strap for the thumb and the Selection® Easy-to-Wrap Closure System. Adjust the resistance by using only one stay for mild resistance and two for increased resistance. The lacing system helps the orthosis to conform around the wrist and give a superior fit.

Made with the airy and comfortable SELECTION material with all the functionality this incorporates.

Recommended Range Of Application

Lateral epicondylitis, hypermobility, rheumatic problems and unspecific wrist pain.

| Item No. | Product Description | Size | Measurement | L/R |
|-----------|---------------------|---------|-------------|-------|
| 352081010 | Selection Bromsa | X-Small | <11 cm | Left |
| 352081011 | Selection Bromsa | Small | 11 - 14 cm | Left |
| 352081012 | Selection Bromsa | Medium | 14 - 17 cm | Left |
| 352081013 | Selection Bromsa | Large | 17 - 20 cm | Left |
| 352081014 | Selection Bromsa | X-Large | 20 - 23 cm | Left |
| 352082010 | Selection Bromsa | X-Small | <11 cm | Left |
| 352082011 | Selection Bromsa | Small | 11 - 14 cm | Right |
| 352082012 | Selection Bromsa | Medium | 14 - 17 cm | Right |
| 352082013 | Selection Bromsa | Large | 17 - 20 cm | Right |
| 352082014 | Selection Bromsa | X-Large | 20 - 23 cm | Right |

