

# SWASH® STEADY

(Replaces previous SWASH Classic. A “Conversion Chart” can be found below under the DOCUMENTS tab.)

Swash is a unique dynamic hip brace that allows controlled variable abduction and dynamic movement in all three planes of motion. The goal of SWASH® is to prevent excessive hip adduction when sitting, standing, and walking as usual in cerebral palsy. However, perhaps the most important benefit the SWASH® offers is guidance for hip alignment.

SWASH Steady offers maximum trunk control, the choice when patient is primarily non-ambulatory (GMFCS IV-V), has low trunk tone and/or very limited trunk control strength or when greater pelvic coverage triggers more desirable neuro-sensory motor response. In sitting, the base is widened for support, which helps the child to get an upright balanced posture and at the same time stretch the adductor muscles.

A complete SWASH orthosis requires a Pelvic Section, Thigh Cuffs, and Uprights. Please order each item separately.

SWASH Steady Pelvic Section has a polyethylene pelvic section that extends from L2/L3 to distal margin of the sacrum, providing maximum thoracic and pelvic support.

Uprights are available in 2 different versions:

115° = narrow sitting base / 123° = wide sitting base

## Recommended Range Of Application

When maximum trunk control is required, i.e., candidates who lack muscle strength or upper body control to sit upright. When patient is primarily non-ambulatory (GMFCS IV - V). When the greater pelvic coverage triggers more desirable neuro-sensory motor response.

## Measurements

At level of the natural waist for maximum trunk stability.



Item No.	Product Description	Size	Circumference	Length	Diameter
288520001	STEADY Pelvic Section	1	40-46.5 cm	-	-
288520002	STEADY Pelvic Section	2	46-55 cm	-	-
288520003	STEADY Pelvic Section	3	54-63 cm	-	-
288520004	STEADY Pelvic Section	4	62-72 cm	-	-
288540601	STEADY Thigh Cuff	1	21-25 cm	-	-
288540802	STEADY Thigh Cuff	2	25-29 cm	-	-
288540803	STEADY Thigh Cuff	3	29-33 cm	-	-
288540804	STEADY Thigh Cuff	4	33-38 cm	-	-
288540805	STEADY Thigh Cuff	5	38-44 cm	-	-
288540806	STEADY Thigh Cuff	6	44-50 cm	-	-
288550601	STEADY Uprights 115°	1	-	31 cm	6 mm
288550802	STEADY Uprights 115°	2	-	37.5 cm	8 mm
288550803	STEADY Uprights 115°	3	-	44 cm	8 mm
288550804	STEADY Uprights 115°	4	-	48 cm	8 mm
288550805	STEADY Uprights 115°	5	-	54 cm	8 mm
288550806	STEADY Uprights 115°	6	-	60 cm	8 mm
288560601	STEADY Uprights 123°	1	-	31 cm	6 mm
288560802	STEADY Uprights 123°	2	-	37.5 cm	8 mm
288560803	STEADY Uprights 123°	3	-	44 cm	8 mm
288560804	STEADY Uprights 123°	4	-	48 cm	8 mm
288560805	STEADY Uprights 123°	5	-	54 cm	8 mm
288560806	STEADY Uprights 123°	6	-	60 cm	8 mm