

NEW

allard|USA

MultiMotion

Low Profile Contracture Management Joints

Now...with Integrated Uprights!!!



CORRECTIVE

Low Profile, Adjustable Torque
Adjustable ROM, Locking Lever



FREE MOTION

Position on opposite side of corrective
joint to provide desired torsional
stiffness to the orthosis

Save Costs with these new **Fabrication Time-Saving** options!

allard|USA

For detailed product and ordering information:
Call **888-678-6548** or Visit www.allardusa.com

The MultiMotion dynamic corrective system is a combination of two components – the corrective joint and the free motion or ROM joint. All components are available in sizes Small and Regular and may be ordered separately. T-Bars are available with choice of 12mm or 22mm offset.



CORRECTIVE JOINTS


FREE MOTION JOINTS




Bar Adapter - Straight or w/T-Bar



Bar Adapter - Straight or w/T-Bar



Integrated Uprights – Straight or w/T-Bar



Integrated Uprights – Straight or w/T-Bar

ROM JOINTS



Quick MultiMotion Joint and Upright Item No. Summary						
MULTIMOTION JOINTS	Standard		W/ T-Bar 12 mm Offset		W/ T-Bar 22 mm Offset	
Size:	Small	Regular	Small	Regular	Small	Regular
CORRECTIVE	69200 0000	69100 0000	69213 0000	69113 0000	69214 0000	69114 0000
CORRECTIVE W/ UPRIGHTS	69122 0000	69121 0000	69125 0000	69123 0000	69126 0000	69124 0000
FREE MOTION	69505 0000	69500 0000	69518 0000	69513 0000	69519 0000	69514 0000
FREE MOTION W/ UPRIGHTS	69128 0000	69127 0000	69131 0000	69129 0000	69132 0000	69130 0000
ROM	69555 0000	69550 0000	N/A	N/A	N/A	N/A
MULTIMOTION UPRIGHTS	Straight		Offset		T-Bar	
Bar Width:	Small	Regular	Small	Regular	Small	Regular
Low Rigidity	69725 0000	69726 0000	69727 0000	69728 0000	69729 0000	69724 0000
High Rigidity	69701 0000	69702 0000	69705 0000	69706 0000	69722 0000	69723 0000



For detailed product and ordering information:
Call **888-678-6548** or Visit **www.allardusa.com**