

Vission™ ROM Wrap Style Knee Support

Our Vission™ ROM Wrap Knee Support combines range of motion control with medial/lateral stability and compressive support for mild to moderate knee instabilities.

Measurements

Measurement: Knee circumference.



Item No.	Product Description	Size	Measurement
435621	Vission™ ROM Wrap Style Knee Support	X-Small	12in - 13in
435631	Vission™ ROM Wrap Style Knee Support	Small	13in - 14in
435641	Vission™ ROM Wrap Style Knee Support	Medium	14in - 15in
435651	Vission™ ROM Wrap Style Knee Support	Large	15in - 17in
435661	Vission™ ROM Wrap Style Knee Support	XLarge	17in - 19in