

Post-Op Knee Brace

Our Post-Op Knee Brace is a post-surgical knee brace that provides stabilization with controlled range of motion during rehabilitation.

Flexion stops: 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°, 105°, 120°
Extension stops: 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°

Recommended Range Of Application

Post-operative fixation and/or range of motion control. Short term stabilization after knee injury.



Item No.	Size	Product description
475150000	One Size	Vission™ Post-Op Knee Brace