SELECTION Wrist Orthosis Rigid



Selection Wrist Orthosis Rigid has the same features as the soft version but includes a thumb section for increased thumb support and stability. The thumb strap can be individually adjusted. It comes with a pre-shaped plastic spoon on the palmar side and a soft narrow plastic spoon on the dorsal side. This product is intended to be used when the first CMC-joint and wrist needs support and off-loading in cases such as rheumatic disorders, soft tissue injuries or attritional wear.

Made with the airy and comfortable SELECTION material with all the functionality this incorporates.

Recommended Range Of Application
Inflammation, work related strain or when rheumatic and/or degenerative changes are present. Also, for soft tissue injuries and attritional wear.

Contraindications
Non known.

Measurements
Wrist circumference.

Item No.	Product Description	Size	Measurement	L/R	Color
352061010	Rigid	X-Small	< 11 cm	Left	Svart
352061011	Rigid	Small	11 - 14 cm	Left	Svart
352061012	Rigid	Medium	14 - 17 cm	Left	Svart
352061013	Rigid	Large	17 - 20 cm	Left	Svart
352061014	Rigid	X-Large	20 - 23 cm	Left	Svart
352061015	Rigid	XX-Large	23 - 26 cm	Left	Svart
352062010	Rigid	X-Small	< 11 cm	Right	Svart
352062011	Rigid	Small	11 - 14 cm	Right	Svart
352062012	Rigid	Medium	14 - 17 cm	Right	Svart
352062013	Rigid	Large	17 - 20 cm	Right	Svart
352062014	Rigid	X-Large	20 - 23 cm	Right	Svart
352062015	Rigid	XX-Large	23 - 26 cm	Right	Svart